

Calibration of Condensation Particle Counter Detection Efficiency Using Mono-Disperse Aerosols

Derek Oberreit, Siqin He

Kanomax FMT



Presentation Outline

- What is a CPC and how is it used for liquid metrology
- CPC Principle of Operation
- Overview of existing calibration methodology
- Description of novel calibration method
- Data discussion


What is a Condensation Particle Counter (CPC)?

- Particles in air are grown to an optically detectable diameter by condensation of a known working vapor
- Condensation happens when:
 - Working vapor is Supersaturated i.e. Vapor pressure, P_{vap} is greater than Saturation Vapor Pressure, P_{sat} at a given gas temperature.
 - The particle size is large enough to overcome surface tension driven evaporation.
- Can be coupled with a nebulizer-aerosolizer to detect particles in liquids

$$Diameter_{Threshold} = \frac{\gamma V_{vap}}{k_B T} \left[\ln \left(\frac{P_{vap}}{P_{sat}} \right) \right]^{-1}$$

$$S = \frac{P_{vap}}{P_{sat}} \propto \frac{1}{D_t}$$


Your video will go here



How does a CPC operate?

- Saturator
 - Aerosol passes through conduit with warm walls wetted with working fluid
- Condenser
 - Saturated aerosol passes laminarily through a conduit with cold walls. For working vapors with Lewis number > 1 , vapor will become supersaturated
 - Peak supersaturation, S_{max} , increases with the difference between Condenser and Saturator temperatures
 - Particles larger than a threshold diameter will grow to an optically detectable size
- Threshold Scanning
 - Condenser is cycled through several set temperatures
 - Each temperature provides a new minimum detected size
 - Temperatures need calibration to ensure consistent operation between devices

Your video will go here



$$Le = \frac{\alpha \text{ (thermal diffusivity, air)}}{D \text{ (mass diffusivity, vapor)}}$$

The CPC Transfer Function

$$N = \int_0^{\infty} \frac{dn}{dD_p} \eta_A(D_p, T_{sat}, T_{cond}) \eta_T(D_p, Q, L) dD_p$$


- N is the concentration measured by the CPC
- n is the true particle concentration
- D_p is the particle diameter
- η_A is the fraction of particles that are activated (i.e. grown to a larger size through condensation)
- η_T is the fraction of particles that are not lost to the conduit walls via diffusion
- T_{sat} and T_{cond} are the saturator and condenser temperatures
- Q is the CPC flow
- L is the effective length of the transport conduit prior to activation (after which transport losses are negligible).

Objective: Determine T_{cond} to achieve a target η_A for a given diameter

Your video will go here

How are CPC traditionally calibrated?


- Aerosol Generation and Conditioning
 - Polydisperse aerosol generated using atomizer, tube furnace, or electro spray
 - Charge conditioners apply a known charge distribution to the aerosol
- Aerosol size selection
 - Differential Mobility Analyzer (DMA) isolates particles of a desired size
 - Only charged particles exit the DMA
- Reference Detector (Electrometer)
 - Faraday Cage Aerosol Electrometer measures the current from the charged particles and converts to a concentration



R. A. Fletcher, G. W. Mulholland, M. R. Winchester, R. L. King & D. B. Klinedinst (2009)
Calibration of a Condensation Particle Counter Using a NIST Traceable Method,
Aerosol Science and Technology, 43:5, 425-441, DOI: 10.1080/02786820802716735

Why are traditional calibration methods tricky?

- Aerosol Generation and Conditioning
 - Polydisperse aerosols skew the challenge aerosol
 - Imperfect charge conditioners skew selection and detection
- Aerosol size selection
 - Transfer function of DMA introduces error due to the range of particle sizes that exit at a given voltage
 - DMAs allow large, multiply charged particles to pass through
- Reference Detector
 - Faraday cage aerosol electrometers unable to distinguish multiply charged particles from singly charged
 - Transport losses in sample tubing and within the reference detector introduce an offset



Goals for a new approach to Calibrating CPCs

- Repeatable
 - Improve agreement between instruments
- Accurate
 - Eliminate artifacts e.g. Charge, Transport, Polydispersity
 - Prove measurement capability
- Practical
 - Limit the complexity of the instrumentation
 - Reduce potential errors

Objective: Determine T_{cond} to achieve a target η_A for a given diameter

The new approach to Calibrating CPCs




- Use monodisperse particles for challenge aerosol
 - Eliminates polydispersity artifacts
 - Calibration diameters are limited to available materials
- Use a “boosted” CPC as reference detector
 - Tuned for 100% detection of all challenge particle sizes
 - Eliminates charge artifacts
 - Allows for lower challenge concentrations compared to electrometers
- Decouple the transport and activation efficiencies
 - Avoids correcting for transport losses
 - Improves repeatability and accuracy

Objective: Determine T_{cond} to achieve a target η_A for a given diameter

Your video will go here

Selecting the challenge aerosol


- Proteins used because they are naturally a fixed molecular weight and are readily available
- Cytochrome-C ~3.3 nm, Immunoglobulin G (IgG) ~ 8nm, Immunoglobulin M (IgM) ~ 15nm



Your video will go here

The Experimental Setup

- Challenge particles are aerosolized
 - Large proteins dissolved in UPW
 - Three solutions prepared and are paired with a multi-channel peristaltic pump
 - Size distribution measured by boosted CPC to verify system operation
- Aerosol passes through an Annular Flow Ion Mobility Classifier (AFIMC)
 - Same principle of operation as a Differential Mobility Analyzer
 - Classifier is used to isolate monomers from multimers
 - Transfer function can be relatively broad
 - Boosted FastCPC used as reference




Your video will go here

Collecting the data

- Software controls pump and classifier voltage
 - Select vial and flowrate
 - Set the classifier voltage to the isolate the particle size being aerosolized
- Software adjusts condenser temperature of CPC under calibration and reads concentration data from both CPCs
 - Temperature range is selected to extend from $\eta_A < 10\%$ to 100% (the asymptote)
- Following a measurement software moves to next vial/size

Detection Efficiency Calibration (Raw)




Analyzing the data

Your video will go here

- Determine the Transport Factor and scale data
 - Ratio of transport efficiencies between test and ref: $TF = \eta_{T,test} / \eta_{T,ref}$
 - $(N_{test} / N_{ref}) / TF = 1$.
- Determine T_{cond} for each size
 - To set the 50% detection efficiency find T_{cond} where $(N_{test} / N_{ref}) / TF = 0.5$.
 - 50% detection selected to minimize bias to smaller sizes. Note: industry standard is to specify minimum detected size (1%-10% detection efficiency)

Detection Efficiency Calibration (Corrected)



Key Takeaways

- Traditional methods for calibrating operating temperatures for Condensation Particle Counters have several sources of error
 - Differences in transport efficiencies between reference and test instrument
 - Skewing due to the shape of the challenge aerosol distribution
 - Over-counting due to multiply charged particles
- Calibrating CPCs using monodisperse particles and temperature stepping mitigates these sources of error

Your video will go here

Thank you to the R&D team at Kanomax FMT for their contributions to many aspects of this work

Derek Oberreit, Kanomax FMT

Derek.Oberreit@kanomaxfmt.com

Note: principle concepts of the described method are patent pending

